


THE PATENTED
PENTADENTATE CHELATOR (PDC)
OF STRUCTURE




is the best chelator
 in Immobilized Metal ion Affinity
 Chromatography (IMAC)
 for
 6xhis-tagged protein purification
 and
 Zn-protein purification
 and
 General purification of proteins


Tran Quang Minh, PhD
 Managing director
 Affiland s.a.
 Liège - Belgium
 October, 2006

INITIAL ATTEMPTS
 at
Purification of a soluble
6xhis-tagged protein
 using
 the patented PDC kit from AFFILAND
 (Cu-PDC, Zn-PDC, Ni-PDC & Co-PDC)
 in comparison with
 the Ni-NTA gel from QIAGEN
 by
 EUROGENTEC S.A. and col.
 DELPHI genetics
 (Belgium)
 June, 2006


M : Markers expressed in Dalton
 BE : Crude extract = lysate
 FL : Flow through
 Buffer E : 100 mM Imidazole pH 7.5
 Buffer F : 300 mM Imidazole pH 7.5




Ni




Ni : Qiagen



Zn



Cu



Co



Immobilized PDC

patented since 1998 by Affiland

Matrix---N(CH₂COOH)-CH₂-CH₂- CH₂-CH₂-CH(COOH)N(CH₂-COOH)₂


See details at <http://www.affiland.com/imac/pdc.htm>

CHARACTERISTICS

- At basic pHs : **Immobilized PDC** is a pentadentate chelator (**Figure 1**).
- At physiological pHs : **Immobilized PDC** is a mixture of tetradentate and pentadentate chelator (**Figure 2 + Figure 1**)

RESULTS

- 1) **Immobilized PDC** can form stable **octahedral complexes M²⁺- PDC** with metal ions **M²⁺** such as Cu²⁺ , Zn²⁺ , Ni²⁺ , Co²⁺ , etc...
- 2) Such **octahedral complexes** possess **only one free site** for interaction with electron donor group of DNA, proteins... (**Figure 1 & Figure 2**).
- 3) As **histidine** and **histidine residue** are known as the ONLY α-amino-acid and its corresponding residue capable of forming stable **octahedral complexes** with **M²⁺**; the **M²⁺- PDC** with **one free site** cited in point 1) and 2) results so in **high selective binding for histidine residues** of 6xhis-tagged proteins and/or histidine containing **proteins**.
- 4) The presence of **only one free site** of the **M²⁺- PDC** complexes for interaction also results in **low concentration of imidazole (50-100 mM)** at physiological pH being sufficient for the elution of 6xhis-tagged proteins from a **M²⁺- PDC** column.
- 5) Due to steric effects of the interaction between a **histidine residue** and the **free site of the M²⁺- PDC** complexes, the **patented PDC kit** composed of 4 columns of Cu-PDC, Zn-PDC, Ni-PDC and Co-PDC is the ONLY WAY of saving time and money during initial attempts at purification of a his-tagged protein.



Solid matrix

Figure 1




Figure 2

Immobilized NTA (Nitrilotriacetic acid) chelator

Ref.: Hochuli et al.; J. Chromatogr. 411, 177 (1987)


Agarose---NH- (CH₂)_n-CH(COOH)N(CH₂-COOH)₂ (n = 2 or 4)

CHARACTERISTICS

- At basic pHs : **Immobilized NTA** is a tetradentate chelator (**Figure 3**).
- At physiological pHs : **Immobilized NTA** is a mixture of tridentate and tetradentate chelator (**Figure 3 + Figure 4**)


RESULTS

- 1) **Immobilized NTA** can form octahedral complexes with metal ions such as Ni²⁺, Co²⁺ ...
- 2) Such complexes leave **two free sites** for interaction with histidine residues of the protein to be purified (**Figure 3 & Figure 4**). The selectivity for histidine residue binding may not be optimal because of the **second free site which can be occupied by other amino-acid residues** of the protein of interest such as cysteine, lysine, tryptophan, **once the first free site is occupied by a histidine residue**.
- 3) The presence of the tridentate NTA complex form (**Figure 4**) at physiological pHs may provoke an important leaching of Ni²⁺ ions during chromatography steps.
- 4) In addition, the **two free sites** cited above results in the need of **high concentration of imidazole (250 mM or more) at physiological pH** for the elution of 6xhis-tagged proteins from a **Ni²⁺ - NTA or Co²⁺ - NTA** column.
- 5) Due to steric effects in interaction between **histidine residue** and **metal chelate NTA complex**, the Ni-NTA and/or Co-NTA CANNOT be the UNIVERSAL WAY for general purification of his-tagged proteins (see enclosed example).
- 6) Contrary to the description given by many manufacturers, the **- NH group of histidine residue is known to be incapable of coordinating with the Ni²⁺ or Co²⁺ ions of the Ni-chelator complexes**.



Solid matrix

Figure 3



Solid matrix

Figure 4

Immobilized IDA (Iminodiacetic acid) chelator

Ref.: Porath et al. - Nature, 258, 589 (1975); Biochem, 22, 1621 (1983)

Sephacrose---N(CH₂-COOH)₂

CHARACTERISTICS

- **Immobilized IDA** is normally a mixture of bidentate and tridentate chelator at physiological pH. It is therefore incapable of forming octahedral complexes with M²⁺ ions. It can **ONLY** form a square or tetrahedral complex with Cu²⁺, Zn²⁺ ... (Figure 5).
- However, **in the presence of electron donor containing cross-linkers, immobilized IDA may form PSEUDO-octahedral complexes with Ni²⁺ ions**, (Figure 6). These PSEUDO-octahedral complexes therefore possess **three free sites** for interaction with the electron donor groups of the proteins to be purified (Figure 5 & Figure 6).

RESULTS

- 1) **NON-selective binding for 6xhis-tagged proteins because the Ni²⁺-IDA is simply a PSEUDO-octahedral complex while Ni²⁺- histidine residue is a TRUE octahedral complex.** The his-tagged proteins purified from Ni²⁺- IDA column cannot be pure and will be contaminated by a lot of other proteins.
- 2) **The three free sites cited above result in the need for a very high concentration of imidazole (500 mM or more) at physiological pH for the elution of 6xhis-tagged proteins.** Indeed, these **CHAOTROPIC CONDITIONS of elution** may result in the activity of such purified proteins being partially or totally destroyed.

Remarks :

- A very high affinity constant is not recommended for Affinity chromatography because the elution step of the product to be purified in this case, will necessitate CHAOTROPIC CONDITIONS i.e. very acid pHs or a very important concentration of organic compounds).
- The **HIGH AFFINITY** of a protein for an affinity matrix **DOES NOT MEAN** the **HIGH CAPACITY** and/or **HIGH SELECTIVITY** of the matrix for this protein.

One variant of immobilized IDA:

Immobilized TED {(N,N,N'-Triscarboxymethyl)ethylenediamine}

Resin---NH(CH₂COOH)-CH₂-CH₂-NH(CH₂-COOH)₂

At physiological pHs, the two NH groups of the TED will be partially charged - the presence of the tridentate TED (3 COO⁻) complex form at physiological pHs may provoke an important leaching of Ni²⁺ ions during chromatography steps.

